


# HANDHELD ECG MONITOR



# HANDHELD ECG MONITOR


MD100A1



MD100B



OEM & Distributor are welcome



- Exquisite, compact & comfortable design
- Easy carry, Simple operation, Quick result
- Large and high resolution LCD screen with back light
- Fast measurement in just 30 seconds
- Displaying & recording heart rate and waveform
- Auto-evaluation of heart health condition
- Automatic power off design
- USB interface for data transfer
- Accuracy and reliability guaranteed

# HANDHELD ECG MONITOR

## Technical Specifications

Standard Specification	MD100A1	MD100B
<b>Power Supply</b>		
	AAA Battery*2	AAA Battery*2
<b>Lead Position</b>		
	Lead I-Hand/Chest Lead; Lead II or Lead III-Between hands and abdomen or left leg	Lead I-Hand/Chest Lead; Lead II or Lead III-Between hands and abdomen or left leg
<b>Display</b>		
	240X160 dot-matrix single color LCD	240X160 dot-matrix single color LCD
<b>Memory</b>		
	30 Seconds x 200 Measured ECG Waveforms	30 Seconds x 200 Measured ECG Waveforms
<b>Waveform Output Method</b>		
	LCD and USB to PC	LCD and USB to PC
<b>ECG Bandwidth</b>		
	0.5Hz-75Hz(+0.4dB, -3dB)	1Hz-75Hz(-3dB)
<b>ADC Resolution</b>		
	10Bit	10Bit
<b>Sampling Rate</b>		
	200Hz	250Hz
<b>Heart Rate Measurement Range</b>		
	30bpm-240bpm	30bpm-240bpm
<b>Accuracy of Heart Rate Measurement</b>		
	30bpm-100bpm: ≤2bpm; 101bpm-240bpm: ±4bpm	30bpm-100bpm: ≤2bpm; 101bpm-240bpm: ±4bpm
<b>Operating Temperature</b>		
	5°C~40°C	5°C~40°C
<b>Storage Temperature</b>		
	-20°C~40°C	-20°C~55°C
<b>Operating Humidity</b>		
	≤80%RH	≤80%RH
<b>Storage Humidity</b>		
	≤93%RH	≤93%RH
<b>Dimensions</b>		
	136MM(L) x84mm(W) x21mm(H)	136MM(L) x84mm(W) x21mm(H)
<b>Other Functions (Optional)</b>	<b>MD100A1</b>	<b>MD100B</b>
<b>Data USB Cable (including CD-ROM)</b>		
	Transferring data records to the computer	Transferring data records to the computer
<b>Data USB Cable (including CD-ROM)</b>		
	For the sake of high quality ECG, cable measurement can be used	For the sake of high quality ECG, cable measurement can be used

AUTHORIZED AGENT



**Innovative Medical Marketing of VA, LLC**

Earl A. Aversano, B.S., M.B.A. --- CEO  
354 McLaws Circle, Suite 1  
Williamsburg, VA 23185  
888-317-1808

info@immva.com  
www.immva.com

**Provide for better life**